

# DIFFERENTIAL PEACK EXPRESIONS USING CHIP SEQUENCING (CHIP-SEQ) ANALYSIS

RIMA ZINJUWADIA

MSC SEM-4

BIOINFORMATICS

# Introduction of ChIP-Seq

- ChIP-seq is a wonderful technique that allows us to interrogate the physical binding interactions between protein and DNA using next-generation sequencing.

# What is chromatin immunoprecipitation?


- Chromatin immunoprecipitation (ChIP) allows us to determine protein-binding sites on DNA. Chromatin is the complex of DNA packaged with histone proteins into nucleosomes. ChIP makes use of reversible cross-links made between DNA and associated proteins by formaldehyde fixation of cells or tissue.
- The fixed chromatin is physically sheared and DNA fragments associated with a particular protein are selectively immunoprecipitated and analysed. Analysis can be on a locus-by-locus basis using PCR, but more commonly ChIP is interrogated with microarrays (ChIP-chip) or next-generation sequencing (ChIP-seq)

# How does ChIP-seq work?

- Chromatin immunoprecipitation sequencing, or ChIP-seq, combines ChIP with next-generation sequencing .
- ChIP-seq protocols have been adapted from ChIP-chip methods: proteins are cross-linked to their bound DNA by formaldehyde treatment, cells are homogenized, and chromatin is sheared and immunoprecipitated with antibody-bound magnetic beads.
- The immunoprecipitated DNA is then used as the input for a next-generation sequencing library prep protocol, where it is sequenced and analysed for DNA binding sites.

# How does ChIP-seq work?

- Chromatin immunoprecipitation sequencing, or ChIP-seq, combines ChIP with next-generation sequencing.
- ChIP-seq protocols have been adapted from ChIP-chip methods: proteins are cross-linked to their bound DNA by formaldehyde treatment, cells are homogenized, and chromatin is sheared and immunoprecipitated with antibody-bound magnetic beads.
- The immunoprecipitated DNA is then used as the input for a next-generation sequencing library prep protocol, where it is sequenced and analysed for DNA binding sites.



# Advantages of ChIP-Seq

- Captures DNA targets for transcription factors or histone modifications across the entire genome of any organism
- Defines transcription factor binding sites
- Reveals gene regulatory networks in combination with RNA sequencing and methylation analysis
- Offers compatibility with various input DNA samples

# Computational analysis of ChIP-Seq

- As with many high-throughput sequencing approaches, ChIP-seq generates extremely large data sets, for which appropriate computational analysis methods are required. To predict DNA-binding sites from ChIP-seq read count data, peak calling methods have been developed. The most popular method is MACS which empirically models the shift size of ChIP-Seq tags, and uses it to improve the spatial resolution of predicted binding sites.



# Selection of ChIP-Seq Data (NCBI-SRA)

## Treated SRR708018

[SRX4010478](#): GSM3120639: 3D\_H\_input; Homo sapiens; ChIP-Seq

1 ILLUMINA (Illumina HiSeq 4000) run: 17.6M spots, 882.3M bases, 326.6Mb downloads

Submitted by: NCBI (GEO)

Study: Effects of culture conditions on epigenomic profiles of brain tumor cells

[PRJNA454151](#) • [SRP143840](#) • [All experiments](#) • [All runs](#)

[hide Abstract](#)

We performed epigenomic analysis of brain tumor cells that were collected from micro-engineered three-dimensional tumor models. We used a low-input epigenomic analysis method known as microfluidic-oscillatory-washing-based chromatin immunoprecipitation with sequencing (MOWChIP-seq) to analyze genome-wide histone modification (H3K4me3). We compared H3K4me3 patterns in standard 2D cultures and 3D cultures based on type I collagen hydrogels, under both normoxic and hypoxic conditions. Our work illustrates a direct connection between cell culture or tissue niche condition and genome-wide alterations in histone modification. Overall design: We obtained genome-wide H3K4me3 profiles in U251 cells cultured under different conditions (3D vs. 2D, Hypoxia vs. Normoxia). The MOWChIP-seq experiments were performed using 1000 cells per assay as described in our previous publication (Cao et al. Nature Methods 12 (2015) 959-962). We generated two replicates (R1 and R2) for each sample.

Sample: 3D\_H\_input

[SAMN08998685](#) • SRS3232109 • [All experiments](#) • [All runs](#)

Organism: [Homo sapiens](#)

Library:

Instrument: Illumina HiSeq 4000

Strategy: ChIP-Seq

Source: GENOMIC

Selection: ChIP

Layout: SINGLE

**Construction protocol:** The day after U251s were seeded in the collagen hydrogels, media was refreshed and well plates or flasks were placed in an incubator at normoxic (21% O<sub>2</sub>) or hypoxic (1% O<sub>2</sub>) conditions. Samples were incubated for 72 hours. Cells in flasks were trypsinized, rinsed once in ice cold PBS, then re-suspended in ice cold PBS at 1,000 cell/μl. Cells seeded in collagen hydrogels were obtained by digesting the collagen in a solution of 0.5% collagenase (Thermo Fisher, Waltham, MA) and 1% FBS in Hanks Buffered Salt Solution (HBSS)( Lonza). Collagen scaffolds were submerged in collagenase solution and incubated at 37°C and 5% CO<sub>2</sub> for 2 hours. Digested collagen solution was collected, rinsed once in ice cold PBS, and re-suspended in ice cold PBS at 1,000 cell/μl. 1 μl of phenylmethyldifonyl fluoride (PMSF) (Sigma-Aldrich, St. Louis, MO) and 1 μl of protease inhibitor cocktail (1X concentration) (Sigma-Aldrich, St. Louis, MO) were added to a 100 μl aliquot of cell suspension for each culture condition. All ChIP-seq libraries were constructed using Accel-NGS 2S plus DNA library kit (Swift Bioscience) following the manufacturer's instructions.

Experiment attributes:

[GEO Accession](#): GSM3120639

Links:

Runs: 1 run, 17.6M spots, 882.3M bases, [326.6Mb](#)

Run	# of Spots	# of Bases	Size	Published
<a href="#">SRR7080719</a>	17,645,919	882.3M	326.6Mb	2019-02-04

## Control SRR7080719

[SRX4010477](#): GSM3120638: 3D\_H\_R2; Homo sapiens; ChIP-Seq

1 ILLUMINA (Illumina HiSeq 4000) run: 11.4M spots, 570M bases, 214.8Mb downloads

Submitted by: NCBI (GEO)

Study: Effects of culture conditions on epigenomic profiles of brain tumor cells

[PRJNA454151](#) • [SRP143840](#) • [All experiments](#) • [All runs](#)

[hide Abstract](#)

We performed epigenomic analysis of brain tumor cells that were collected from micro-engineered three-dimensional tumor models. We used a low-input epigenomic analysis method known as microfluidic-oscillatory-washing-based chromatin immunoprecipitation with sequencing (MOWChIP-seq) to analyze genome-wide histone modification (H3K4me3). We compared H3K4me3 patterns in standard 2D cultures and 3D cultures based on type I collagen hydrogels, under both normoxic and hypoxic conditions. Our work illustrates a direct connection between cell culture or tissue niche condition and genome-wide alterations in histone modification. Overall design: We obtained genome-wide H3K4me3 profiles in U251 cells cultured under different conditions (3D vs. 2D, Hypoxia vs. Normoxia). The MOWChIP-seq experiments were performed using 1000 cells per assay as described in our previous publication (Cao et al. Nature Methods 12 (2015) 959-962). We generated two replicates (R1 and R2) for each sample.

Sample: 3D\_H\_R2

[SAMN08998686](#) • SRS3232108 • [All experiments](#) • [All runs](#)

Organism: [Homo sapiens](#)

Library:

Instrument: Illumina HiSeq 4000

Strategy: ChIP-Seq

Source: GENOMIC

Selection: ChIP

Layout: SINGLE

**Construction protocol:** The day after U251s were seeded in the collagen hydrogels, media was refreshed and well plates or flasks were placed in an incubator at normoxic (21% O<sub>2</sub>) or hypoxic (1% O<sub>2</sub>) conditions. Samples were incubated for 72 hours. Cells in flasks were trypsinized, rinsed once in ice cold PBS, then re-suspended in ice cold PBS at 1,000 cell/μl. Cells seeded in collagen hydrogels were obtained by digesting the collagen in a solution of 0.5% collagenase (Thermo Fisher, Waltham, MA) and 1% FBS in Hanks Buffered Salt Solution (HBSS)( Lonza). Collagen scaffolds were submerged in collagenase solution and incubated at 37°C and 5% CO<sub>2</sub> for 2 hours. Digested collagen solution was collected, rinsed once in ice cold PBS, and re-suspended in ice cold PBS at 1,000 cell/μl. 1 μl of phenylmethyldifonyl fluoride (PMSF) (Sigma-Aldrich, St. Louis, MO) and 1 μl of protease inhibitor cocktail (1X concentration) (Sigma-Aldrich, St. Louis, MO) were added to a 100 μl aliquot of cell suspension for each culture condition. All ChIP-seq libraries were constructed using Accel-NGS 2S plus DNA library kit (Swift Bioscience) following the manufacturer's instructions.

Experiment attributes:

[GEO Accession](#): GSM3120638

Links:

Runs: 1 run, 11.4M spots, 570M bases, [214.8Mb](#)

Run	# of Spots	# of Bases	Size	Published
<a href="#">SRR7080718</a>	11,399,827	570M	214.8Mb	2019-02-04




# Single End

- In single-end reading, the sequencer reads a fragment from only one end to the other, generating the sequence of base pairs. In this less sequencing is required and it is used for general purposes like differential expression analysis. Single-read sequencing can be a good choice for certain methods such as small RNA-Seq or chromatin immunoprecipitation sequencing (ChIP-Seq).

Single-end reads



# Workflow of the ChIP-Seq Analysis



# FastQC

- Modern high throughput sequencers can generate hundreds of millions of sequences in a single run. Before analysing this sequence to draw biological conclusions you should always perform some simple quality control checks to ensure that the raw data looks good and there are no problems or biases in your data which may affect how you can usefully use it.
- Most sequencers will generate a QC report as part of their analysis pipeline, but this is usually only focused on identifying problems which were generated by the sequencer itself. FastQC aims to provide a QC report which can spot problems which originate either in the sequencer or in the starting library material.
- FastQC can be run in one of two modes. It can either run as a stand alone interactive application for the immediate analysis of small numbers of FastQ files, or it can be run in a non-interactive mode where it would be suitable for integrating into a larger analysis pipeline for the systematic processing of large numbers of files.

# FastQC result of Control - SRR7080719

## FastQC Report


### Summary

- ✓ Basic Statistics
- ✓ Per base sequence quality
- ✓ Per sequence quality scores
- ⚠ Per base sequence content
- ⚠ Per sequence GC content
- ✓ Per base N content
- ✓ Sequence Length Distribution
- ✓ Sequence Duplication Levels
- ✓ Overrepresented sequences
- ✓ Adapter Content
- ✗ Kmer Content


### Basic Statistics

Measure	Value
Filename	Control_SRR7080719.fastq
File type	Conventional base calls
Encoding	Sanger / Illumina 1.9
Total Sequences	17645919
Sequences flagged as poor quality	0
Sequence length	50
%GC	41


### Per base sequence quality




### ✓ Per sequence quality scores




### ① Per sequence GC content




### ① Per base sequence content




### ✓ Per base N content



## Sequence Length Distribution




## Adapter Content




## Overrepresented sequences

No overrepresented sequences

## Sequence Duplication Levels



## Kmer Content





# FastQC result of Treated - SRR7080718

## FastQC Report


### Summary

- ✓ [Basic Statistics](#)
- ✓ [Per base sequence quality](#)
- ✓ [Per sequence quality scores](#)
- ⚠ [Per base sequence content](#)
- ⚠ [Per sequence GC content](#)
- ✓ [Per base N content](#)
- ✓ [Sequence Length Distribution](#)
- ✓ [Sequence Duplication Levels](#)
- ✓ [Overrepresented sequences](#)
- ✓ [Adapter Content](#)
- ✗ [Kmer Content](#)

### ✓ Basic Statistics


Measure	Value
Filename	treated_SRR7080718.fastq
File type	Conventional base calls
Encoding	Sanger / Illumina 1.9
Total Sequences	11399827
Sequences flagged as poor quality	0
Sequence length	50
%GC	39

### ✓ Per base sequence quality







### ✓ Per sequence quality scores




### 1 Per base sequence content




### 1 Per sequence GC content




### ✓ Per base N content



## Sequence Length Distribution




## Sequence Duplication Levels




## Overrepresented sequences

No overrepresented sequences

## Adapter Content



## Kmer Content



# MAPPING ON REFERENCE GENOME

---

Use Bowtie2 to align reads.

<http://bowtie-bio.sourceforge.net/bowtie2/manual.shtml#options>

# Bowtie 2

- Bowtie 2 is an ultrafast and memory-efficient tool for aligning sequencing reads to long reference sequences. It is particularly good at aligning reads of about 50 up to 100s of characters to relatively long (e.g. mammalian) genomes.
- Bowtie 2 outputs alignments in SAM format, enabling interoperability with a large number of other tools (e.g. SAMtools, GATK) that use SAM. Bowtie 2 is distributed under the GPLv3 license, and it runs on the command line under Windows, Mac OS X and Linux

# Bowtie 2

- Bowtie 2 is often the first step in pipelines for comparative genomics, including for variation calling, ChIP-seq, RNA-seq, BS-seq.
- Bowtie uses indexed genome for the alignment in order to keep its memory footprint small. Because of time constraints we will build the index only for one chromosome of the human genome. For this we need the chromosome sequence in fasta format.

# Bowtie 2

## Steps:

- Indexing using Bowtie2
- Mapping using Bowtie2
- Align the Control\_SRR7080719 reads using Bowtie2
- Align the treated\_SRR7080718 reads using Bowtie2
- Converting SAM to BAM using samtools

# BAM File Formate

- A BAM file (\*.bam) is the compressed binary version of a SAM file that is used to represent aligned sequences up to 128 Mb. SAM and BAM formats are described in detail at <https://samtools.github.io/hts-specs/SAMv1.pdf>.
- BAM files use the file naming format of SampleName\_S#.bam, where # is the sample number determined by the order that samples are listed for the run. In multi-node mode, the S# is set to S1, regardless the order of the sample.
- BAM files contain a header section and an alignment section.
- Header —Contains information about the entire file, such as sample name, sample length, and alignment method. Alignments in the alignments section are associated with specific information in the header section.
- Alignments—Contains read name, read sequence, read quality, alignment information, and custom tags. The read name includes the chromosome, start coordinate, alignment quality, and the match descriptor string



# SAM file format

- SAM stands for Sequence Alignment/Map format. It is a TAB-delimited text format consisting of a header section, which is optional, and an alignment section. If present, the header must be prior to the alignments. Header lines start with '@', while alignment lines do not. Each alignment line has 11 mandatory fields for essential alignment information such as mapping position, and variable number of optional fields for flexible or aligner specific information.

# The alignment section: mandatory fields

Col	Field	Type	Regexp/Range	Brief description
1	QNAME	String	[!-?A-~]{1,254}	Query template NAME
2	FLAG	Int	$[0, 2^{16} - 1]$	bitwise FLAG
3	RNAME	String	\* [:rname:^*=] [:rname:]*	Reference sequence NAME <sup>9</sup>
4	POS	Int	$[0, 2^{31} - 1]$	1-based leftmost mapping POSition
5	MAPQ	Int	$[0, 2^8 - 1]$	MAPping Quality
6	CIGAR	String	\* ([0-9]+[MIDNSHPX=])+	CIGAR string
7	RNEXT	String	\* = [:rname:^*=] [:rname:]*	Reference name of the mate/next read
8	PNEXT	Int	$[0, 2^{31} - 1]$	Position of the mate/next read
9	TLEN	Int	$[-2^{31} + 1, 2^{31} - 1]$	observed Template LENgth
10	SEQ	String	\* [A-Za-z=.]+	segment SEQUENCE
11	QUAL	String	[!-~]+	ASCII of Phred-scaled base QUALity+33

# The alignment section: optional fields

Type	Regex matching VALUE	Description
A	[!~]	Printable character
i	[-+]?[0-9]+	Signed integer <sup>12</sup>
f	[-+]?[0-9]*\.[0-9]+([eE] [-+]?[0-9]+)?	Single-precision floating number
Z	[!~]*	Printable string, including space
H	([0-9A-F] [0-9A-F])*	Byte array in the Hex format <sup>13</sup>
B	[cCsSiIf] (, [-+]?[0-9]*\.[0-9]+([eE] [-+]?[0-9]+)?)*	Integer or numeric array

# SAMtools

- Samtools is a set of utilities that manipulate alignments in the BAM format. It imports from and exports to the SAM (Sequence Alignment/Map) format, does sorting, merging and indexing, and allows to retrieve reads in any regions swiftly.
- Samtools is designed to work on a stream. It regards an input file '-' as the standard input (stdin) and an output file '-' as the standard output (stdout). Several commands can thus be combined with Unix pipes. Samtools always output warning and error messages to the standard error output (stderr).
- Samtools is also able to open a BAM (not SAM) file on a remote FTP or HTTP server if the BAM file name starts with 'ftp://' or 'http://'. Samtools checks the current working directory for the index file and will download the index upon absence. Samtools does not retrieve the entire alignment file unless it is asked to do so.

# SAMtools

## Steps:

- sort sam file into bam sorted
- Remove duplicate reads
- Only pick reads which is mapping uniquely
- Indexing of .bam

# PEAK CALLING

---

MACS<sub>2</sub>

# MACS<sub>2</sub>

- With the improvement of sequencing techniques, chromatin immunoprecipitation followed by high throughput sequencing (ChIP-Seq) is getting popular to study genome-wide protein-DNA interactions.
- To address the lack of powerful ChIP-Seq analysis method, we present a novel algorithm, named Model-based Analysis of ChIP-Seq (MACS), for identifying transcript factor binding sites.
- MACS captures the influence of genome complexity to evaluate the significance of enriched ChIP regions, and MACS improves the spatial resolution of binding sites through combining the information of both sequencing tag position and orientation. MACS can be easily used for ChIP-Seq data alone, or with control sample with the increase of specificity.



# MOTIF-ANALYSIS

---

MEME-ChIPtool

# MEME-ChIP


- In genetics, a sequence motif is a nucleotide or amino-acid sequence pattern that is widespread and has, or is conjectured to have, a biological significance. For proteins, a sequence motif is distinguished from a structural motif, a motif formed by the three-dimensional arrangement of amino acids which may not be adjacent.
- 'Motif discovery' (or 'motif finding') in biological sequences can be defined as the problem of finding short similar sequence elements (building the 'motif') shared by a set of nucleotide or protein sequences with a common biological function.

# MEME-ChIP

- All the peaks fasta sequences retrieved for the motif analysis can be used online for MEME-Chip tool for further motif analysis.
- MEME-Chip performs motif discovery, motif enrichment analysis and clustering on large nucleotide datasets.
- First go to MEME suite from your browser

# The MEME Suite

Motif-based sequence analysis tools



**MEME**  
Multiple EM for Motif Elicitation

**CentriMo**  
Local Motif Enrichment Analysis

**FIMO**  
Find Individual Motif Occurrences

**DREME**  
Discriminative Regular Expression Motif Elicitation

**AME**  
Analysis of Motif Enrichment

**MAST**  
Motif Alignment & Search Tool

**MEME-ChIP**  
Motif Analysis of Large Nucleotide Datasets

**SpaMo**  
Species Motif Analysis Tool

**MCAST**  
Motif Cluster Alignment and Search Tool

**GLAM2**  
Gapped Local Alignment of Motifs

**GOMo**  
Gene Ontology for Motifs

**GLAM2Scan**  
Scanning with Gapped Motifs

**Tomtom**  
Motif Comparison Tool

**GT-Scan**  
Identifying Unique Genomic Targets

**MoMo**  
Modification Motifs

MEME Suite 4.12.0

**Motif Discovery**

MEME

DREME

MEME-ChIP

GLAM2

MoMo

**Motif Enrichment**

CentriMo

AME

SpaMo

GOMo

**Motif Scanning**

FIMO

MAST

MCAST

GLAM2Scan

**Motif Comparison**

Tomtom

**Manual**

OVERVIEW

Motif Discovery

MEME

DREME

MEME-ChIP

GLAM2

MoMo

Motif Enrichment

Go to Motif discovery and select for MEME-CHIP and upload the peak fasta sequences file at the section "Input the primary sequences" and click on "start search"

← → ↻ ① meme-suite.org/tools/meme-chip

**MEME Suite 4.12.0**

▼ Motif Discovery

MEME  
DREME  
MEME-CHIP  
GLAM2  
MoMo

► Motif Enrichment

► Motif Scanning

► Motif Comparison

► Manual

► Guides & Tutorials

► Sample Outputs

► File Format Reference

► Databases

► Download & Install

► Help

► Alternate Servers

► Authors & Citing

▼ Recent Jobs


Thursday 31 August

MEME-CHIP 17:45 ✕

MEME 16:07 ✕

MEME 16:01 ✕

DREME 15:50 ✕



**MEME-CHIP**  
Motif Analysis of Large Nucleotide Datasets  
Version 4.12.0

Data Submission Form

Perform motif discovery, motif enrichment analysis and clustering on large nucleotide datasets.

**Select the motif discovery and enrichment mode**

☒ Normal mode ☐ Discriminative mode ?

**Select the sequence alphabet**

Use sequences with a standard alphabet or specify a custom alphabet. ?

☒ DNA, RNA or Protein ☐ Custom  No file chosen

**Input the primary sequences**

Enter the (equal-length) nucleotide sequences to be analyzed. ?

No file chosen ?

**Input the motifs**

Select, upload or enter a set of known motifs. ?

?

?

**Input job details**

(Optional) Enter your email address. ?

(Optional) Enter a job description. ?






► Universal options

# Results from MEME-CHIP :

## MOTIFS

The significant motifs (E-value  $\leq 0.05$ ) found by the programs MEME, DREME and CentriMo; clustered by similarity and ordered by E-value.

Expand All Clusters Collapse All Clusters

Motif Found	Discovery/Enrichment Program	E-value	Known or Similar Motifs	Distribution	SpaMo & FIMO
 <p>Reverse Complement ⇌</p>	<a href="#">DREME</a>	2.1e-454		Not Centrally Enriched	<ul style="list-style-type: none"> <li><a href="#">Motif Sites in GFF3</a></li> </ul>
 <p>Reverse Complement ⇌ Show 1 More ⓘ</p>	<a href="#">DREME</a>	1.1e-364		Not Centrally Enriched	<ul style="list-style-type: none"> <li><a href="#">Motif Sites in GFF3</a></li> </ul>
 <p>Reverse Complement ⇌ Show 1 More ⓘ</p>	<a href="#">DREME</a>	2.4e-272		Not Centrally Enriched	<ul style="list-style-type: none"> <li><a href="#">Motif Sites in GFF3</a></li> </ul>
 <p>Reverse Complement ⇌</p>	<a href="#">MEME</a>	2.9e-091			<ul style="list-style-type: none"> <li><a href="#">Motif Sites in GFF3</a></li> </ul>



Reverse Complement ⇄ Show 5 More ⓘ CentriMo Group ↻ ⓘ

[CentriMo](#)

4.3e-006

[POU3F4 \(MA0789.1\)](#)



- [Motif Sites in GFF3](#)

Motif Found

Discovery/Enrichment Program ⓘ

E-value ⓘ

Known or Similar Motifs ⓘ

Distribution ⓘ

SpaMo & FIMO ⓘ



Reverse Complement ⇄ Show 1 More ⓘ CentriMo Group ↻ ⓘ

[CentriMo](#)

1.6e-004

[Meox2\\_DBD](#)



- [Motif Sites in GFF3](#)

Motif Found

Discovery/Enrichment Program ⓘ

E-value ⓘ

Known or Similar Motifs ⓘ

Distribution ⓘ

SpaMo & FIMO ⓘ



Reverse Complement ⇄

[CentriMo](#)

6.8e-003

[Arid3a \(MA0151.1\)](#)



- [Motif Sites in GFF3](#)

Motif Found

Discovery/Enrichment Program ⓘ

E-value ⓘ

Known or Similar Motifs ⓘ

Distribution ⓘ

SpaMo & FIMO ⓘ



Reverse Complement ⇄

[CentriMo](#)

2.3e-002

[POU2F3\\_DBD\\_2](#)



- [Motif Sites in GFF3](#)

Motif Found

Discovery/Enrichment Program ⓘ

E-value ⓘ

Known or Similar Motifs ⓘ

Distribution ⓘ

SpaMo & FIMO ⓘ



Reverse Complement ⇄

[CentriMo](#)

4.6e-002

[Hoxd9\\_DBD\\_3](#)



- [Motif Sites in GFF3](#)



# CHIP PEAK ANNOTATION, COMPARISON, AND VISUALIZATION

ChIPseeker

# Chip Peak Annotation, Comparison, And Visualization

- Chromatin immunoprecipitation followed by high-throughput sequencing (ChIP-seq) has become standard technologies for genome wide identification of DNA-binding protein target sites. After read mappings and peak callings, the peak should be annotated to answer the biological questions.
- Annotation also create the possibility of integrating expression profile data to predict gene expression regulation. ChIPseeker was developed for annotating nearest genes and genomic features to peaks.
- ChIP peak data set comparison is also very important. ChIPseeker (Yu, Wang, and He 2015) support statistical testing of significant overlap among ChIP seq data sets, and incorporate open access database GEO for users to compare their own dataset to those deposited in database. Converting genome coordinations from one genome version to another is also supported, making this comparison available for different genome version and different species.


# Chip Peak Annotation, Comparison, And Visualization

- Several visualization functions are implemented to visualize the coverage of the ChIP seq data, peak annotation, average profile and heatmap of peaks binding to TSS region.
- Functional enrichment analysis of the peaks can be performed by my Bioconductor packages DOSE(Yu et al. 2015), ReactomePA(Yu and He 2016), clusterProfiler(Yu et al. 2012).


# ChIPseeker

- ChIPseeker is an R package for annotating ChIP-seq data analysis. It supports annotating ChIP peaks and provides functions to visualize ChIP peaks coverage over chromosomes and profiles of peaks binding to TSS regions. Comparison of ChIP peak profiles and annotation are also supported. Moreover, it supports evaluating significant overlap among ChIP-seq datasets.


# Results generated after running ChIPseeker tool:




covplot




peakHeatmap




**plotAnnoPie**



**plotAvgProf2**



**tagHeatmap**



**upsetplot**